


## Product datasheet for **GA104340**

### Human Syndecan 4 (SDC4) activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | SDC4 CRISPRa kit - CRISPR gene activation of human syndecan 4           |
| Symbol:       | SDC4  |
| UniProt ID:   | <a href="#">P31431</a>  |
| Synonyms:     | SYND4   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

#### Product images:



[View online >](#)