

Product datasheet for **GA104350**

Human SEC14L1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SEC14L1 CRISPRa kit - CRISPR gene activation of human SEC14 like lipid binding 1
Symbol:	SEC14L1
Components:	GA104350G1, SEC14L1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA104350G2, SEC14L1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA104350G3, SEC14L1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q92503
Synonyms:	PRELID4A; SEC14L
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)