

Product datasheet for **GA104353**

Human SEL1L activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SEL1L CRISPRa kit - CRISPR gene activation of human SEL1L adaptor subunit of ERAD E3 ubiquitin ligase
Symbol:	SEL1L
UniProt ID:	Q9UBV2
Synonyms:	PRO1063; SEL1-LIKE; SEL1L1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)