

Product datasheet for **GA104361**

Human SEPP1 (SELENOP) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SELENOP CRISPRa kit - CRISPR gene activation of human selenoprotein P
Symbol:	SELENOP
UniProt ID:	<u>P49908</u>
Synonyms:	SELP; SeP; SEPP; SEPP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)