

Product datasheet for **GA104371**

Human SF2 (SRSF1) activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | SRSF1 CRISPRa kit - CRISPR gene activation of human serine and arginine rich splicing factor 1 |
| Symbol: | SRSF1 |
| UniProt ID: | Q07955 |
| Synonyms: | ASF; SF2; SF2p33; SFRS1; SRp30a |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)