

## Product datasheet for **GA104373**

### Human SFRS3 (SRSF3) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	SRSF3 CRISPRa kit - CRISPR gene activation of human serine and arginine rich splicing factor 3
Symbol:	SRSF3
UniProt ID:	<a href="#">P84103</a>
Synonyms:	SFRS3; SRp20
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)