

Product datasheet for **GA104374**

Human SFRS4 (SRSF4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SRSF4 CRISPRa kit - CRISPR gene activation of human serine and arginine rich splicing factor 4
Symbol:	SRSF4
UniProt ID:	Q08170
Synonyms:	SFRS4; SRP75
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)