

Product datasheet for **GA104392**

Human SH3BGRL activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | SH3BGRL CRISPRa kit - CRISPR gene activation of human SH3 domain binding glutamate rich protein like |
| Symbol: | SH3BGRL |
| UniProt ID: | <u>O75368</u> |
| Synonyms: | HEL-S-115; SH3BGR |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)