

Product datasheet for **GA104401**

Human SHFM3 (FBXW4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	FBXW4 CRISPRa kit - CRISPR gene activation of human F-box and WD repeat domain containing 4
Symbol:	FBXW4
Components:	GA104401G1, SHFM3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA104401G2, SHFM3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA104401G3, SHFM3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P57775</u>
Synonyms:	DAC; FBW4; FBWD4; SHFM3; SHSF3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)