


Product datasheet for **GA104403**

Human SHMT1 activation kit by CRISPRa

Product data:

| | |
|---------------|---------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits |
| Description: | SHMT1 CRISPRa kit - CRISPR gene activation of human serine hydroxymethyltransferase 1 |
| Symbol: | SHMT1 |
| UniProt ID: | P34896 |
| Synonyms: | CSHMT; SHMT |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)