

Product datasheet for **GA104408**

Human SIAH1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SIAH1 CRISPRa kit - CRISPR gene activation of human siah E3 ubiquitin protein ligase 1
Symbol:	SIAH1
UniProt ID:	<u>Q8IUQ4</u>
Synonyms:	SIAH1A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)