

## Product datasheet for **GA104409**

### Human SIAH2 activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | SIAH2 CRISPRa kit - CRISPR gene activation of human siah E3 ubiquitin protein ligase 2 |
| Symbol:       | SIAH2  |
| UniProt ID:   | <u><a href="#">O43255</a></u>  |
| Synonyms:     | hSiah2   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                |

#### Product images:



[View online >](#)