

Product datasheet for **GA104416**

Human Melanoma gp100 (PMEL) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PMEL CRISPRa kit - CRISPR gene activation of human premelanosome protein
Symbol:	PMEL
UniProt ID:	P40967
Synonyms:	D12S53E; gp100; ME20; ME20-M; ME20M; P1; P100; PMEL17; SI; SIL; SILV
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)