

Product datasheet for **GA104433**

Human SLC4A3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SLC4A3 CRISPRa kit - CRISPR gene activation of human solute carrier family 4 member 3
Symbol:	SLC4A3
Components:	GA104433G1, SLC4A3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA104433G2, SLC4A3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA104433G3, SLC4A3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P48751</u>
Synonyms:	AE3; CAE3/BAE3; SLC2C
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)