

Product datasheet for **GA104445**

Human Band 3 (SLC4A1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SLC4A1 CRISPRa kit - CRISPR gene activation of human solute carrier family 4 member 1 (Diego blood group)
Symbol:	SLC4A1
UniProt ID:	P02730
Synonyms:	AE1; BND3; CD233; CHC; DI; EMPB3; EPB3; FR; RTA1A; SAO; SPH4; SW; WD; WD1; WR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)