

Product datasheet for **GA104449**

Human Smoothelin (SMTN) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SMTN CRISPRa kit - CRISPR gene activation of human smoothelin
Symbol:	SMTN
Components:	GA104449G1, Smoothelin gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA104449G2, Smoothelin gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA104449G3, Smoothelin gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P53814</u>
Synonyms:	FLJ35168; FLJ35365; FLJ35620; FLJ38597
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)