

Product datasheet for **GA104469**

Human **SLC8A3** activation kit by **CRISPRa**

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | SLC8A3 CRISPRa kit - CRISPR gene activation of human solute carrier family 8 member A3 |
| Symbol: | SLC8A3 |
| UniProt ID: | P57103 |
| Synonyms: | NCX3 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)