

Product datasheet for **GA104519**

Human SMN1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SMN1 CRISPRa kit - CRISPR gene activation of human survival of motor neuron 1, telomeric
Symbol:	SMN1
UniProt ID:	Q16637
Synonyms:	BCD541; GEMIN1; SMA; SMA1; SMA2; SMA3; SMA4; SMA@; SMN; SMNT; T-BCD541; TDRD16A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)