


Product datasheet for **GA104532**

Human SNAPC3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SNAPC3 CRISPRa kit - CRISPR gene activation of human small nuclear RNA activating complex polypeptide 3
Symbol:	SNAPC3
Components:	GA104532G1, SNAPC3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA104532G2, SNAPC3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA104532G3, SNAPC3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q92966</u>
Synonyms:	PTFbeta; SNAP50
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)