

Product datasheet for **GA104546**

Human SNRPD3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SNRPD3 CRISPRa kit - CRISPR gene activation of human small nuclear ribonucleoprotein D3 polypeptide
Symbol:	SNRPD3
UniProt ID:	P62318
Synonyms:	Sm-D3; SMD3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)