

Product datasheet for **GA104550**

Human SNRPN activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SNRPN CRISPRa kit - CRISPR gene activation of human small nuclear ribonucleoprotein polypeptide N
Symbol:	SNRPN
UniProt ID:	P63162
Synonyms:	HCERN3; PWCR; PWS; RT-LI; SM-D; sm-N; SMN; SNRNP-N; SNURF-SNRPN
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)