

Product datasheet for **GA104559**

Human Superoxide Dismutase 3 (SOD3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SOD3 CRISPRa kit - CRISPR gene activation of human superoxide dismutase 3
Symbol:	SOD3
UniProt ID:	P08294
Synonyms:	EC-SOD
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)