

Product datasheet for **GA104579**

Human Transcription factor Sp4 (SP4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SP4 CRISPRa kit - CRISPR gene activation of human Sp4 transcription factor
Symbol:	SP4
Components:	GA104579G1, Transcription factor Sp4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA104579G2, Transcription factor Sp4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA104579G3, Transcription factor Sp4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q02446</u>
Synonyms:	HF1B; SPR-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)