


Product datasheet for **GA104581**

Human SPAG1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SPAG1 CRISPRa kit - CRISPR gene activation of human sperm associated antigen 1
Symbol:	SPAG1
Components:	GA104581G1, SPAG1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA104581G2, SPAG1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA104581G3, SPAG1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q07617</u>
Synonyms:	CILD28; CT140; HEL-S-268; HSD-3.8; SP75; TPIS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)