

Product datasheet for **GA104587**

Human SPG7 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | SPG7 CRISPRa kit - CRISPR gene activation of human SPG7 matrix AAA peptidase subunit, paraplegin |
| Symbol: | SPG7 |
| UniProt ID: | Q9UQ90 |
| Synonyms: | CAR; CMAR; PGN; SPG5C |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)