

## Product datasheet for **GA104599**

### Human **SPRR1B** activation kit by **CRISPRa**

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | SPRR1B CRISPRa kit - CRISPR gene activation of human small proline rich protein 1B                                                                                                                                                                                                                            |
| Symbol:       | SPRR1B                                                                                                                                                                                                                                                                                                        |
| Components:   | GA104599G1, SPRR1B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA104599G2, SPRR1B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA104599G3, SPRR1B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">P22528</a></u>                                                                                                                                                                                                                                                                                 |
| Synonyms:     | CORNIFIN; GADD33; SPRR1                                                                                                                                                                                                                                                                                       |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)