

Product datasheet for **GA104639**

Human **SSR2** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SSR2 CRISPRa kit - CRISPR gene activation of human signal sequence receptor subunit 2
Symbol:	SSR2
UniProt ID:	P43308
Synonyms:	HSD25; TLAP; TRAP-BETA; TRAPB
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)