


Product datasheet for **GA104640**

Human SSR3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SSR3 CRISPRa kit - CRISPR gene activation of human signal sequence receptor subunit 3
Symbol:	SSR3
Components:	GA104640G1, SSR3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA104640G2, SSR3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA104640G3, SSR3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9UNL2</u>
Synonyms:	TRAPG
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)