

Product datasheet for **GA104646**

Human Somatostatin Receptor 3 (SSTR3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SSTR3 CRISPRa kit - CRISPR gene activation of human somatostatin receptor 3
Symbol:	SSTR3
UniProt ID:	P32745
Synonyms:	SS-3-R; SS3-R; SS3R; SSR-28
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)