

## Product datasheet for **GA104697**

### Human Suppressor of Ty 4 homolog 1 (SUPT4H1) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	SUPT4H1 CRISPRa kit - CRISPR gene activation of human SPT4 homolog, DSIF elongation factor subunit
Symbol:	SUPT4H1
UniProt ID:	<a href="#">P63272</a>
Synonyms:	SPT4; SPT4H; Supt4a; SUPT4H
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)