

Product datasheet for **GA104700**

Human SUPV3L1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SUPV3L1 CRISPRa kit - CRISPR gene activation of human Suv3 like RNA helicase
Symbol:	SUPV3L1
Components:	GA104700G1, SUPV3L1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA104700G2, SUPV3L1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA104700G3, SUPV3L1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8IYB8</u>
Synonyms:	SUV3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)