

Product datasheet for **GA104724**

Human Neurokinin A Receptor (TACR2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TACR2 CRISPRa kit - CRISPR gene activation of human tachykinin receptor 2
Symbol:	TACR2
UniProt ID:	P21452
Synonyms:	NK2R; NKNAR; SKR; TAC2R
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)