

Product datasheet for **GA104725**

Human Neurokinin B (TAC3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TAC3 CRISPRa kit - CRISPR gene activation of human tachykinin 3
Symbol:	TAC3
UniProt ID:	<u>Q9UHF0</u>
Synonyms:	gamma tachykinin 3; neurokinin-B protein; neurokinin B-like protein; neurokinin beta; neuromedin K, neurokinin beta; NKB; NKB, NKNB, PRO1155, ZNEUROK1; NKNB; preprotachykinin B; PRO1155; tachykinin 3; ZNEUROK1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)