

Product datasheet for **GA104728**

Human Neurokinin 1 Receptor (TACR1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TACR1 CRISPRa kit - CRISPR gene activation of human tachykinin receptor 1
Symbol:	TACR1
UniProt ID:	P25103
Synonyms:	NK1R; NK1R; SPR; TAC1R
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)