

Product datasheet for **GA104729**

Human Neurokinin 3 Receptor NK3 (TACR3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TACR3 CRISPRa kit - CRISPR gene activation of human tachykinin receptor 3
Symbol:	TACR3
UniProt ID:	P29371
Synonyms:	HH11; NK-3R; NK3R; NKR; TAC3RL
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)