

Product datasheet for **GA104732**

Human TAF2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TAF2 CRISPRa kit - CRISPR gene activation of human TATA-box binding protein associated factor 2
Symbol:	TAF2
UniProt ID:	Q6P1X5
Synonyms:	CIF150; MRT40; TAF2B; TAFII150
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)