

## Product datasheet for **GA104733**

### Human TAF4 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TAF4 CRISPRa kit - CRISPR gene activation of human TATA-box binding protein associated factor 4
Symbol:	TAF4
UniProt ID:	<a href="#">O00268</a>
Synonyms:	TAF2C; TAF2C1; TAF4A; TAFII130; TAFII135
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:

