

Product datasheet for **GA104734**

Human TAF4B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TAF4B CRISPRa kit - CRISPR gene activation of human TATA-box binding protein associated factor 4b
Symbol:	TAF4B
UniProt ID:	Q92750
Synonyms:	SPGF13; TAF2C2; TAFII105
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)