

Product datasheet for **GA104735**

Human SM22 alpha (TAGLN) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TAGLN CRISPRa kit - CRISPR gene activation of human transgelin
Symbol:	TAGLN
UniProt ID:	Q01995
Synonyms:	SM22; SMCC; TAGLN1; WS3-10
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)