

Product datasheet for **GA104739**

Human TAF9 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TAF9 CRISPRa kit - CRISPR gene activation of human TATA-box binding protein associated factor 9
Symbol:	TAF9
UniProt ID:	Q16594
Synonyms:	MGC:5067; STAF31/32; TAF2G; TAFII-31; TAFII-32; TAFII31; TAFII32; TAFIID32
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)