

Product datasheet for **GA104740**

Human TAF10 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TAF10 CRISPRa kit - CRISPR gene activation of human TATA-box binding protein associated factor 10
Symbol:	TAF10
UniProt ID:	Q12962
Synonyms:	TAF2A; TAF2H; TAFII30
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)