

Product datasheet for **GA104741**

Human TAF11 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TAF11 CRISPRa kit - CRISPR gene activation of human TATA-box binding protein associated factor 11
Symbol:	TAF11
UniProt ID:	Q15544
Synonyms:	MGC:15243; PRO2134; TAF2I; TAFII28
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)