

Product datasheet for **GA104743**

Human TAF13 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TAF13 CRISPRa kit - CRISPR gene activation of human TATA-box binding protein associated factor 13
Symbol:	TAF13
UniProt ID:	Q15543
Synonyms:	TAF(II)18; TAF2K; TAFII-18; TAFII18
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)