

Product datasheet for **GA104753**

Human TARS1 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | TARS CRISPRa kit - CRISPR gene activation of human threonyl-tRNA synthetase |
| Symbol: | TARS1 |
| UniProt ID: | P26639 |
| Synonyms: | ThrRS |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)