

Product datasheet for **GA104782**

Human TCF3 / E2A (TCF3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TCF3 CRISPRa kit - CRISPR gene activation of human transcription factor 3
Symbol:	TCF3
UniProt ID:	P15923
Synonyms:	AGM8; bHLHb21; E2A; E47; ITF1; TCF-3; VDIR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)