

Product datasheet for **GA104789**

Human ZNF354A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF354A CRISPRa kit – CRISPR gene activation of human zinc finger protein 354A
Symbol:	ZNF354A
Components:	GA104789G1, ZNF354A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA104789G2, ZNF354A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA104789G3, ZNF354A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>O60765</u>
Synonyms:	EZNF; HEL104; HKL1; KID-1; KID1; TCF17
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

