

Product datasheet for **GA104799**

Human TCP11 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TCP11 CRISPRa kit - CRISPR gene activation of human t-complex 11
Symbol:	TCP11
Components:	GA104799G1, TCP11 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA104799G2, TCP11 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA104799G3, TCP11 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8WWU5
Synonyms:	D6S230E; FPPR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

