

## Product datasheet for **GA104806**

### Human DYNLT3 activation kit by CRISPRa

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                               |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| <b>Description:</b>  | DYNLT3 CRISPRa kit - CRISPR gene activation of human dynein light chain Tctex-type 3                                                                                                                                                                                                                          |
| <b>Symbol:</b>       | DYNLT3                                                                                                                                                                                                                                                                                                        |
| <b>Components:</b>   | GA104806G1, DYNLT3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA104806G2, DYNLT3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA104806G3, DYNLT3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| <b>UniProt ID:</b>   | <u>P51808</u>                                                                                                                                                                                                                                                                                                 |
| <b>Synonyms:</b>     | RP3; TCTEIL; TCTEXIL                                                                                                                                                                                                                                                                                          |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:

