

Product datasheet for **GA104814**

Human TEF1 (TEAD1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TEAD1 CRISPRa kit - CRISPR gene activation of human TEA domain transcription factor 1
Symbol:	TEAD1
UniProt ID:	P28347
Synonyms:	AA; NTEF-1; REF1; TCF-13; TCF13; TEAD-1; TEF-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)