

Product datasheet for **GA104815**

Human TEAD4 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | TEAD4 CRISPRa kit - CRISPR gene activation of human TEA domain transcription factor 4 |
| Symbol: | TEAD4 |
| UniProt ID: | Q15561 |
| Synonyms: | EFTR-2; hrTEF-1B; RTEF1; TCF13L1; TEF-3; TEF3; TEFR-1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)